


What is the Brain Imaging Data Structure and why you should know about this!


07/11/2024 Julius Welzel







What is the Brain Imaging Data Structure and why you should know about this!

The Past, Present, and Future of BIDS



The Past, Present, and Future of BIDS





What is BIDS?



About 20+ different EEG systems record in different file formats



Aims of BIDS




→ Makes (neuroimaging-)data FAIR

DATA STRUCTURE

- BIDS is based on simple file formats and folder structures
- Provides a "easy-to-understand" folder structure
- Provides a naming convention for files and directories
- Specifies file formats to use for modalities

METADATA

- Prevents metadata getting lost
- Some metadata is better than no metadata (80/20 rule)



GitHub repository

The Brain Imaging Data Structure


The Brain Imaging Data Structure (BIDS) is a simple and intuitive way to organize and describe data.

A Q Search

This document defines the BIDS specification, which provides many details to help implement the standard. It includes the core specification as well as many extensions to specific brain imaging modalities, and increasingly also to other kinds of data.

If BIDS is new to you, and you would like to learn more about how to adapt your own datasets to match the BIDS specification, we recommend exploring the BIDS Starter Kit. Alternatively, to get started please read the introduction to the specification.

For an overview of the BIDS ecosystem, visit the BIDS homepage. The entire specification can also be downloaded as PDF.



Remi Gau

Remi-Gau

Follow


I wanted to understand the brain, now I tell people how to name files, folders and variables.

My favourite two BIDS principles

Requirement levels

Level of requirement for files and within files:

- 1. REQUIRED
- 2. RECOMMENDED
- 3. OPTIONAL



participant_id	age	group
sub-001	34	С
sub-002	12	Р
sub-003	33	M

BIDS definitions

Dataset - A set of neuroimaging and behavioral data acquired for a purpose of a particular study. A dataset consists of data acquired from one or more subjects, possibly from multiple sessions.

Subject - A person or animal participating in the study. Used interchangeably with term **Participant**.

Session - A logical grouping of neuroimaging and behavioral data consistent across subjects.

Task - A set of structured activities performed by the participant.

Overview files:

- participant.tsv
- scans.tsv
- ...

Subject, session and task specific files

- sub-001_ses-01_task-GoNoGo.eeg
- sub-001_task-Oddball_channels.tsv

Directories and file naming structure

- Data for each subject are organized in subdirectories labeled "sub-<label>", with "<label>" representing the unique identification label of each subject
- If subjects participated in multiple sessions,
 subdirectories labeled "ses-<label>" contain session-specific data within the subject directory
- Each session subdirectory (or the subject subdirectory if no session data exists) contains subdirectories for different data types, only defined if files are present for that type
- For a data file that was collected in a given session from a given subject, the file name MUST begin with the string sub-<label>[_ses-<label>]_task-<label>

```
README.md
dataset_description.json
participants.json
participants.tsv
sub-001
ses-01
eeg/
sub-001_ses-01_task-GoNoGo.eeg
```





Modality agnostic data

dataset_description.json

```
README.md
dataset description.json
     participants.json
     participants.tsv
     sub-001
        ses-01
        ses-02
     sub-002
```

```
{
"Name": "The mother of all experiments",
"BIDSVersion": "1.4.0",
"DatasetType": "raw",
"License": "CCO",
"Authors": [ "Paul Broca", "Carl Wernicke"],
"Funding": [ "NIH F37823MFH1"],
"EthicsApprovals": ["HRPO (Protocol AR0928"],
"DatasetDOI": "10.0.2.3/dfjj.10"
}
```

participants

```
README.md
☐ dataset_description.json
☐ participants.tsv
☐ participants.json
☐ sub-001
☐ ses-01
☐ ses-02
☐ sub-002
```

participants.tsv


```
participant_id age group
sub-001 34 C
sub-002 12 P
sub-003 33 M
```

participants.json

```
"age": {
    "Description": "age of the
    participant", "Units": "years"
    },

"group": {
    "Description": "assigned group",
    "Levels": {
        "C": "control",
        "M": "medication",
        "P": "placebo",
    }
}
```





Modality specific data

Files | eeg

```
____ sub-001
____ ses-01
____ eeg/
```

```
"TaskName": "TASKNAME",
    "SamplingFrequency": 1000,
    "SoftwareFilters": "n/a",
    "EEGChannelCount": 4,
    "EOGChannelCount": 1,
    "EEGReference": "placed on Cz",
    "PowerLineFrequency": 50
}
```

```
eeg/
sub-001_ses-01_task-<label>_eeg.<extension>
sub-001_ses-01_task-<label>_eeg.json
sub-001_ses-01_task-<label>_channels.tsv
sub-001_ses-01_task-<label>_electrodes.tsv
sub-001_ses-01_task-<label>_electrodes.tsv
```

```
impedance
name
                          Z
CP5
        -0.77
                 -0.30
                          0.57
FC5
        -0.77
                 0.30
                          0.57
                                   12
FC1
        -0.29
                 0.31
                          0.91
                                   2
C3
        -0.59
                 0.00
                          0.81
VEOG
        n/a
                 n/a
                          n/a
                                   n/a
```

```
units
                         status
                                  status description
name
        type
CP5
        EEG
                 microV
                                  n/a
                         good
FC5
        EEG
                                  high freq noise
                 microV
                         bad
FC1
        EEG
                 microV
                         good
                                  n/a
C3
        EEG
                 microV
                                  n/a
                         good
VEOG
        EOG
                 microV
                                  n/a
                         good
```

```
{
    "EEGCoordinateSystem": "T1w",
    "EEGCoordinateUnits": "mm",
    "AnatomicalLandmarkCoordinates": {
        "LPA": [-0.067, 1.736e-09, -3.844e-09],
        "NAS": [-4.11e-09, 0.091, -4.541e-10],
        "RPA": [0.064, -6.435e-09, -4.566e-09]
    },
    "AnatomicalLandmarkCoordinateSystem": "T1w",
    "AnatomicalLandmarkCoordinateUnits": "mm",
    "IntendedFor": "sub-01_T1w.nii.gz"
}
```

Metadata | eeg.json

REQUIRED

TaskName, SamplingFrequency,
PowerlineFrequency, EEGReference,
SoftwareFilters

RECOMMENDED

TaskDescription, Manufacturer, EEGChannelCount, EOGChannelCount, ...

OPTIONAL

_

```
TaskName: Seeing stuff,
SamplingFrequency: 2400,
Manufacturer: Brain Products,
ManufacturersModelName: BrainAmp DC,
EEGChannelCount: 87,
EOGChannelCount: 2,
TriggerChannelCount: 1,
PowerLineFrequency: 50,
EEGPlacementScheme: 10 percent system,
EEGReference: single electrode placed on FCz,
SoftwareFilters: {
       Anti-aliasing filter: {
                          half-amplitude cutoff (Hz):
                           500,
                          Roll-off: 6dB/Octave
```

Metadata channels.tsv

REQUIRED

name, type, units

RECOMMENDED

_

OPTIONAL


description, sampling_frequency

name	type	units	description	reference	status	status_description
VEOG	VEOG	υV	left eye	VEOG-, VEOG+	good	n/a
FDI	EMG	υV	left first dorsal interosseous	FDI-, FDI+	good	n/a
Cz	EEG	υV	n/a	REF	bad	high frequency noise
UADC001	MISC	n/a	envelope of audio signal	n/a	good	n/a

Metadata electrodes.tsv

name	X	У	Z	type	material	color
VEOG+	n/a	n/a	n/a	cup	Ag/AgCl	blue
VEOG-	n/a	n/a	n/a	cup	Ag/AgCl	white
FDI+	n/a	n/a	n/a	cup	Ag/AgCl	red
FDI-	n/a	n/a	n/a	cup	Ag/AgCl	red
GND	-0.07	0.00	-0.070	clip-on	Ag/AgCl	pink
Cz	0.00	0.07	0.06	cup	Ag/AgCl	yellow
REF	-0.07	-0.02	-0.01	cup	Ag/AgCl	grey

Column name	Requirement Level	Data type	Description
name	REQUIRED	string	Name of the electrode contact point.
			Values in name MUST be unique.
			This column must appear first in the file.
х	REQUIRED	number	Recorded position along the x-axis.
			This column must appear second in the file.
у	REQUIRED	number	Recorded position along the y-axis.
			This column must appear third in the file.
z	REQUIRED	number or	Recorded position along the z-axis.
		n/a	This column must appear fourth in the file.
type	RECOMMENDED	string	Type of the electrode (for example, cup, ring, clip-on, wire, needle).
			This column may appear anywhere in the file.
material	RECOMMENDED	string	Material of the electrode (for example, Tin, Ag/AgCl, Gold).
			This column may appear anywhere in the file.
impedance	RECOMMENDED	number	Impedance of the electrode, units MUST be in k0hm.
			This column may appear anywhere in the file.
Additional Columns	OPTIONAL	n/a	Additional columns are allowed if they are defined in the associated metadata file.



eeg/

Modality specific data

Files | motion

```
0,2634511
                                     0,092295
                                                0,0086682
                                                           0,9305117
                                                                      0,690106
                          0,694520
                                                                      0,885496
                                     0,1918243
                                                0,8437273
                                                           0,3975710
sub-001/
                          0,0766395
                                     0,2587211
                                                0,5434792
                                                           0,2822837
                                                                      0,2789791
                          0,5779993
                                     0,0456141
                                                0,0490745
                                                           0,9408899
                                                                      0,1533421
    ses-01/
                          0.054556
                                     0,7915927
                                                0,5871733
                                                           0,4669577
                                                                      0,9754468
                          0.966024
                                     0,1962834
                                                0,7114406
                                                           0,3389448
                                                                      0,7194495
       eeg/
                          0,984172
                                     0,5079461
                                                0,1180168
                                                                      0,1753768
                                                           0,7966978
                          0,9883907
                                     0,1557346
                                                0,8002013
                                                           0,6334882
                                                                      0,7526906
  •••
       -motion ∕
          sub-001 ses-01 task-<|abel> tracksys-<|abel> motion.tsv
          sub-001 ses-01 task-<label> tracksys-<label> motion.json
     — sub-001 ses-01 task-<label> channels.tsv
```

```
{
"SamplingFrequency": 60,
"SamplingFrequencyEffective": 60.19,
"TaskName": "BIDS Motion fictive example",
"TrackingSystemName": "imu1",
"TaskDescription": "walking and talking",
"MotionChannelCount": 6,
"SubjectArtefactDescription": "n/a",
"TrackedPointsCount": 2,
"ACCELChannelCount": 3,
"GYROChannelCount": 3,
"Manufacturer": "BWSensing",
"ManufacturersModelName": "BW-imu600",
}
```

0.8098815

0,8952724

0,2326254

0,6683652

0,0480541

0,4384892

0,4886533

0,8529441

```
tracked point
                                                                                                 placement
name
                            component
                                          type
                                                                                   units
                                                                                   m/s^2
                                                                                                 right foot
imul rf acc x
                                         ACCEL
                           Х
                                         ACCEL
                                                                                   m/s^2
                                                                                                 right _foot
imu1_rf_acc_y
                                                       rf
imu1_rf_acc_z
                                         ACCEL
                                                                                   m/s^2
                                                                                                 right _foot
                                                       rf
                           Z
imu1_rf_gyro_x
                                         GYRO
                                                                                   rad/s
                                                                                                 right foot
                           Х
                                                       rf
imu1_rf_gyro_y
                                         GYRO
                                                       rf
                                                                                   rad/s
                                                                                                 right foot
imu1_rf_gyro_z
                                         GYRO
                                                       rf
                                                                                   rad/s
                                                                                                 right _foot
                           Z
```

Metadata motion.json

REQUIRED

TaskName, SamplingFrequency

RECOMMENDED

TaskDescription, RecordingDuration,
MotionChannelCount, <type>ChannelCount,
SoftwareFilters, ...

OPTIONAL

Manufacuturer, RecordingSoftware, ...

```
"SamplingFrequency": 60,
"TaskName": "BIDS Motion fictive example",
"TrackingSystemName": "IMU Right Hand",
"TaskDescription": "walking and talking",
"InstitutionAddress": "Fictive address",
"InstitutionName": "Fictive Institution",
"MotionChannelCount": 18,
"Recording Duration": 4667.641106,
"SubjectArtefactDescription": "n/a",
"TrackedPointsCount": 2,
"ACCELChannelCount": 6,
"GYROChannelCount": 6,
"MAGNChannelCount": 6,
"Manufacturer": "BWSensing",
"ManufacturersModelName": "BW-IMU600",
```

Metadata channels.tsv

REQUIRED

name, component, type, tracked_point, units

RECOMMENDED

placement, reference_frame

OPTIONAL

. . .


status, status_description, sampling_frequency

Restricted keyword list for channel type					
Restricted keyword list for column type in alphabetic order. Note that upper-case is REQUIRED:					
Keyword Description					
Accelerometer channel, one channel for each spatial axis. Column component for the axis MUST be added to the *_channels.tsv file (x, y, or z).					
Angular acceleration channel, one channel for each spatial axis. Column component for the axis MUST be added to the *_channels.tsv file (x, y, or z).					
Gyrometer channel, one channel for each spatial axis. Column component for the axis MUST be added to the *_channels.tsv file (x, y, or z).					
Joint angle channel between two fixed axis belonging to two bodyparts. Angle SHOULD be defined between proximal and distal bodypart in deg.					

name	component	type	tracked_point	units	reference_frame
t1_acc_x	X	ACCEL	LeftFoot	$m/s^{\Lambda}2$	global
t1_acc_y	У	ACCEL	LeftFoot	$m/s^{\Lambda}2$	global
t1_acc_z	Z	ACCEL	LeftFoot	$m/s^{\Lambda}2$	global
t1_gyro_x	X	GYRO	LeftFoot	rad/s	global
t1_gyro_y	У	GYRO	LeftFoot	rad/s	global
t1_gyro_z	Z	GYRO	LeftFoot	rad/s	global

Overview | motion

```
0,26345511
            0,092292015 0,008668652
                                         0,930514317
                                                        0,690193606
                                                                        0,809881135
0,694520294
            0,191824943 0,843726573
                                         0,397571025
                                                        0,88542996
                                                                        0,895276224
0,076633595
            0,258720111 0,547534792
                                         0,282283781
                                                        0,27890791
                                                                        0,232620594
            0,045616941 0,04903375
                                         0,940889749
                                                        0,153318421
                                                                        0,668360752
            0,791513927 0,587116733
                                         0,466957774
                                                        0,975446368
                                                                        0,048053341
            0.196283834 0.711044406
                                         0,338944328
                                                        0,719445195
                                                                        0,438488392
0,98417512
            0,507944361 0,1180168
                                         0,796692478
                                                        0,175376468
                                                                        0,488659533
0,98839607
            0,155737146 0,800206213
                                         0,633481382
                                                        0,752698206
                                                                        0,852943441
```




	{	
		"SamplingFrequency": 60,
		"SamplingFrequencyEffective": 60.19,
		"TaskName": "BIDS Motion fictive example",
		"TrackingSystemName": "imu1",
		"TaskDescription": "walking and talking",
		"MotionChannelCount": 6,
		"RecordingDuration": 10.05,
		"SubjectArtefactDescription": "n/a",
		"TrackedPointsCount" : 2,
		"ACCELChannelCount": 3,
		"GYROChannelCount": 3,
		"Manufacturer": "BWSensing",
		"ManufacturersModelName": "BW-imu600",
	}	
ı		

onset	duration	trial_type
2.45	0.1	heel_strike
3.81	0.1	heel_strike
4.95	0.1	heel_strike
6.11	0.1	heel_strike
7.24	0.1	heel_strike

name	type	units	srate	tracked_point	component	placement
imu1_rf_acc_x	ACCEL	m/s^2	222	rf	х	right_foot
imu1_rf_acc_y	ACCEL	m/s^2	222	rf	У	right _foot
imu1_rf_acc_z	ACCEL	m/s^2	222	rf	Z	right _foot
imu1_rf_gyro_x	GYRO	rad/s	222	rf	Х	right _foot
imu1_rf_gyro_y	GYRO	rad/s	222	rf	У	right _foot
imu1_rf_gyro_z	GYRO	rad/s	222	rf	Z	right _foot





Modality specific data

Task events

```
sub-<label>/
     <data_type>/
          <matches>_events.tsv
          <matches>_events.json
```

```
{
  "channel": {
      "Description": "Channel(s) associated with the event"
},
  "annots": {
      "LongName": "Annotations",
      "Description": "Annotations associated with channels indicated in the channel column.",
      "Levels": {
            "musc": "Muscle artifact. A very common, high frequency, sharp artifact that corresponds with agitation/nervousness in a patient."
      }
    }
}
```

onset	duration	trial_type	response_time	stim_file	channel	annots
1.23	0.65	start	1.435	red_square.jpg	n/a	n/a
5.65	0.65	stop	1.739	blue_square.jpg	n/a	n/a
12.1	2.35	n/a	n/a	n/a	Cz	musc

scans.tsv

```
sub-<label>/
  [ses-<label>/]
  sub-<label>[_ses-<label>]_scans.tsv
```

Column name	Requirement Level	Data type	Description
filename	REQUIRED	string	Relative paths to files. There MUST be exactly one row for each file.
			Values in filename MUST be unique.
			This column must appear first in the file.
acq_time	OPTIONAL	string	Acquisition time refers to when the first data point in each run was acquired. Furthermore, if this header is provided, the acquisition times of all files from the same recording MUST be identical. Datetime format and their anonymization are described in Units.
			This column may appear anywhere in the file.
Additional Columns	OPTIONAL	n/a	Additional columns are allowed.

filename

eeg/sub-control01_task-nback.eeg
eeg/sub-control01_task-walking.eeg
motion/sub-control01_task-walking.tsv

acq_time

1877-06-15T13:45:30

1877-06-15T13:55:33

1877-06-15T13:56:27

*nback.eeg

*walking.eeg

*walking.tsv

What is the Brain Imaging Data Structure and why you should know about this!

For research


Validation and optimizing pipelines

Building a sustainable research environment

Article | Open access | Published: 09 February 2023 **EEG is better left alone**Arnaud Delorme ☑

Scientific Reports 13, Article number

Preprocessing Choices for P3 Analyses with Mobile EEG:A Systematic Literature Review and Interactive Exploration



For you

Load public data from the BIDS format

 EEGLab, FieldTrip and MNE provide importers



Converting your data to BIDS

- Reuse data from your Lab and others
- Foster collaborations



```
sub = \{ '01', '02', '03', '04', '05', '06', '07', '08', '09', '10' \};
% for subject 3 the age is unknown, for subject 2 the sex is not specified
age = [11 96 nan 77 82 87 18 40 26 80];
sex = {'f' [] 'f' 'f' 'm' 'm' 'm' 'm' 'm'}:
for subindx=1:numel(sub)
  cfq = [];
 cfg.datatype = 'eeg';
  % Load vour data
 cfg.dataset = sub{subindx} + '.edf';
  % specify the output directory
 cfg.bidsroot = 'bids';
                = sub{subindx};
  cfg.sub
  % specify the information for the participants.tsv file
 % this is optional, you can also pass other pieces of info
 cfg.participants.age = age(subindx);
  cfg.participants.sex = sex{subindx};
 % provide the mnemonic and long description of the task
                      = 'changedetection':
  cfa.TaskName
  cfg.TaskDescription = 'Subjects were responding as fast as possible upon a
                          change in a visually presented stimulus.';
 % these are EEG specific
 cfg.eeg.PowerLineFrequency = 50; % since recorded in the EU
 cfg.eeg.EEGReference
                            = 'M1'; % left mastoid
  data2bids(cfq):
```

end

BIDS Starterpack

Website

[https://bids-specification.readthedocs.io/en/stable/]

Validator

[https://bids-standard.github.io/bids-validator/]

BIDS Validator v1.14.6

Select a RIDS dataset to validate



Brain Imaging Data Structure v1.9.0







Brain Imaging Data Structure

The BIDS Specification

Introduction

Common principles

Modality agnostic files

Modality specific files

Derivatives

Longitudinal and multi-site studies

Glossary

BIDS Extension Proposals

Appendix

Changelog

The BIDS Starter Kit

Website

Tutorials

GitHub repository

The Brain Imaging Data Structure

The Brain Imaging Data Structure (BIDS) is a simple and intuitive way to organize and describe

This document defines the BIDS specification, which provides many details to help implement the standard. It includes the core specification as well as many extensions to specific brain imaging modalities, and increasingly also to other kinds of data.

If BIDS is new to you, and you would like to learn more about how to adapt your own datasets to match the BIDS specification, we recommend exploring the BIDS Starter Kit. Alternatively, to get started please read the introduction to the specification.

For an overview of the BIDS ecosystem, visit the BIDS homepage. The entire specification can also be downloaded as PDF.







Conclusions








Outlook: EEGManySteps





Thank you for listening carefully



Thanks to **Sein Jeung** for pushing this to completion over the past three years Thanks to all of the **BIDS Maintainers and Devs** who help to implement this Thanks to my **working group**, who have given me the freedom to work in this Thanks to all **participants** who ever provided data which is now in BIDS